

Quality Improvement Guidelines for Bone Tumour Management

A. Gangi · G. Tsoumakidou · X. Buy ·
E. Quoix

Received: 5 October 2009 / Accepted: 5 October 2009 / Published online: 12 February 2010
© Springer Science+Business Media, LLC and the Cardiovascular and Interventional Radiological Society of Europe (CIRSE) 2010

Introduction

Different image-guided percutaneous techniques can be used for treatment or pain palliation in patients with primary or secondary bone tumours. Curative ablation can be applied for the treatment of specific benign or in selected cases of malignant localized bone tumours [1]. Pain palliation therapy of primary and secondary bone tumours [2, 3] can be achieved with safe, fast, effective, and tolerable percutaneous methods [4]. Ablation (chemical, thermal, mechanical), cavitation (radiofrequency ionization), and consolidation (cementoplasty) techniques can be used separately or in combination. Each technique has its indications as well as its own advantages and drawbacks.

Definitions

- *Curative treatment* The therapeutic goal of curative treatment is to completely and definitively ablate the tumour.
- *Palliative treatment* The therapeutic goal of palliative treatment is not to completely ablate the tumour but to achieve pain palliation.

Consolidation Techniques

- *Cementoplasty* Percutaneous injection of methylmethacrylate provides bone strengthening and pain relief

[5, 6]. Cementoplasty is a symptomatic and not a curative treatment [7, 8].

Ablation Techniques

- *Chemical ablation* Injection of ethanol is used to ablate osteolytic bone tumours as well as achieve pain management (neurolysis) [9].
- *Thermal ablation procedures* consist of the following:
 1. *Laser photocoagulation* Thermoablation is produced with near-infrared wavelength lasers (neodymium yttrium aluminum garnet [Nd:YAG] diode laser 800–1000 nm).
 2. *Radiofrequency ablation* Thermoablation is produced by altering the electric current at the tip of the electrode, thus causing local ionic agitation and subsequent frictional heat.
 3. *Cryoablation* The application of extreme cold aims to destroy cells by causing both direct cellular and vascular injury [10, 11].
- *Radiofrequency ionisation (decompression ablation technique)* consists of the following:
 1. Low-temperature bipolar technique: This produces a plasma field at the tip of the electrode that breaks the intermolecular bonds, thus creating a cavity inside the tissue, which leads in turn to decompression.

Treatment Strategy

The management of patients with bone tumours requires consideration of the following:

A. Gangi (✉) · G. Tsoumakidou · X. Buy · E. Quoix
Non-Vascular Interventional Radiology Department, Strasbourg
University Hospital, 67091 Strasbourg, Cedex, France
e-mail: gangi@rad6.u-strasbg.fr

- histology of the tumour with differentiation of benign and malignant tumours;
- careful evaluation of the patient's general condition;
- understanding of the disease process;
- appreciation of the degree of bone destruction (consolidation);
- thorough working knowledge of available treatment options and multidisciplinary decision making; and
- precision of the treatment goal, i.e., curative or palliative.

Tumour Biopsy

Percutaneous image-guided biopsy of the tumour should be obtained when needed. Specific benign bone tumours have characteristic imaging findings (e.g., osteoid osteoma) and thus usually do not require a pretreatment biopsy specimen. For the rest, i.e., benign and primary malignant neoplasms, histologic confirmation is necessary. For chondroblastoma, a biopsy should always be performed to differentiate it from clear cell chondrosarcoma.

Secondary metastatic tumours usually do not require a biopsy specimen. Biopsy should be performed only in cases of unknown primary tumour or multiple primary neoplasms or when there is doubt if the lesion is a secondary localization of the known primary neoplasm.

Curative Percutaneous Bone Tumour Management

Indications

- Benign lesions
 1. Osteoid osteoma;
 2. osteoblastoma (when the maximum diameter is <3cm); and
 3. chondroblastoma.
- Malignant lesions (primary or secondary in nonsurgical patients) are made by multidisciplinary decision. Treatment is reserved for slow-growing cancers with the following typical features:
 1. with <3 proven locations, and
 2. measuring <3 cm in diameter.

Palliative Percutaneous Bone Tumour Management

Indications

In the palliative treatment of painful bone tumours (primary or secondary), the therapeutic goal is not complete ablation of the tumour but rather:

- pain palliation;
- tumour reduction;
- prevention of impending pathologic fractures; and/or
- decompression and debulking of bone tumours (especially for tumours bulging into the spinal canal).

Procedure Set-Up and Sedation

General Considerations

Coagulation studies (prothrombin time, activated partial prothrombin time, international normalized ratio) as well as platelet and white blood counts should be normal. All infections located remotely from the surgical site should be identified and resolved before the intervention is performed. Informed written consent should be obtained before the procedure.

Procedure Set-Up

Dealing with bone tumours requires strict sterile conditions. Antibiotics can be administered intravenously on the day of the procedure to prevent infection.

Sedation

Bone tumour management is usually painful and requires sedation and analgesia. Tumour ablation usually needs deeper sedation, such a general anaesthesia, nerve root blocks, or spinal anaesthesia. For single-level cementoplasty, general anaesthesia is not mandatory; conscious sedation is usually sufficient.

Types of Image Guidance

- *Fluoroscopy* is the most commonly used imaging modality for skeletal interventions, especially when real-time imaging is needed (e.g., cement injection during cementoplasty).
- *Ultrasound* can be used in superficial lytic tumours, i.e., those with a soft-tissue component.
- *Computed axial tomography (CAT) guidance* provides high-resolution images that are useful for tumour ablation (three-dimensional [3-D] imaging is necessary) and help decrease the complication rate.
- *Magnetic resonance imaging (MRI) guidance* is promising for its thermal-monitoring capacity. It also offers the possibility of 3-D imaging, which is necessary for tumour ablation.

- *Flat-panel cone-beam CAT* is becoming more accurate, faster, and therefore very promising.
- *Multimodality imaging*, which is associated with fluoroscopy and CAT scan, allows faster, safer, and more accurate procedures and is particularly useful for complex cases (e.g., spine).

Review of Techniques

Cementoplasty

Patient Selection and Technique

Osteolytic tumours (e.g., metastasis, multiple myeloma, lymphoma) that are located near the vertebral body, acetabulum, and condyles and that cause local pain, disability, and a risk of compression fracture are excellent indications for cementoplasty.

Cementoplasty can be performed with the patient under sedation or general anaesthesia (especially for multilevel vertebroplasty) [12, 13] using fluoroscopic or/and CAT guidance. For vertebroplasty, a 10- or 15-gauge beveled needle (for thoracic/lumbar and cervical level, respectively) should be used. Cement injection should be performed under real-time imaging (fluoroscopic-lateral view, multislice fluoro-CAT).

Patient Outcome

Significant pain relief (>50% decrease in narcotic analgesics used or >50% pain decrease per visual analogue scale) has been reported in 75% to 85% of patients [8, 14–16]. Improved results have been obtained by combining thermal ablation and cementoplasty in cases of paravertebral tumor extension or soft-tissue extension [17].

Complications

Cement leakage has been reported to occur in 38% to 72.5% of patients [18]; lung embolism, infection, and fracture in <1% of patients; and allergic reaction, puncture-site bleeding, and acute increased pain or new areas of pain in 14% of patients [16].

Chemical Ablation with Ethanol

Patient Selection and Technique

Palliative treatment can be offered to patients with osteolytic painful bone tumours without risk of fracture [5]. The procedure should be performed strictly with the patient

under deep sedation or general anaesthesia (due to the pain induced during ethanol injection in bone and soft tissue) and under CAT guidance. According to the size and number of lesions and level of pain relief desired, alcoholization using 3 to 25 ml sterile 95% ethanol is administered either as a multisession technique or as a one-shot technique through a 22-gauge needle.

Nonionic iodinated contrast medium should be always injected into the lesion before ethanol injection to predict the distribution. In patients with large tumours, alcohol should be selectively instilled into regions considered to be responsible for pain (usually the periphery of the tumour and osteolytic areas). The size and shape of the necrosis induced with ethanol is not always reproducible and depends on the degree of vascularization, necrosis, and tissue consistency.

Patient Outcome

Pain relief has been reported in 74% of patients and usually occurs within 24 hours. In 26% of cases, decreased tumour size has been observed, whereas in 18.5% of cases tumour size continues to increase. Duration of pain relief ranges from 10 to 27 weeks. The best results are obtained in patients with small lytic metastases (i.e., diameter ranging from 3 to 6 cm) [9].

Complications

Minor complications include low-grade fever (17%), hyperuricemia, and pain during the first 6 hours. Major complications include infection and accidental leakage causing damage to the contiguous neurologic structures (<1.5%).

Laser Ablation

Patient Selection and Technique

Because of the small size of the ablation zone produced, laser is mostly used for patients with small tumours or in those with contraindications to radiofrequency ablation (RFA) (e.g., metallic implants). Multiple fibres are required for patients with larger tumours.

Osteoid Osteomas

The best indication is tumour size (<1 cm in diameter). The procedure should always be performed with the patient under general anaesthesia or regional block (severe pain is typically experienced during penetration of the nidus). The laser fibre (400–600 µm) is always inserted coaxially into the tumour under CAT guidance through a spinal needle or

drill needle when necessary. The choice of access route (all neurovascular structures should be avoided) and diffusion of heat should be taken into account to avoid complications. If required, additional insulation techniques and thermal monitoring should be applied to protect adjacent vulnerable structures. The positioning of the fibre exactly in the centre of the nidus is the key to a successful procedure. Two-watt power is applied for 6–10 min depending on tumour size. At the end of ablation, 5–10 ml ropivacaine 2 mg/ml is injected (strictly extravascularly) in contact with the periosteum to decrease postprocedural pain.

Painful Bone Metastases

Painful bone metastases can be treated by inserting ≤ 8 simultaneously energized bare fibres into the tumour with approximately 1.5-cm spacing.

Patient Outcome

The success rate for ablation of osteoid osteomas is similar to that of RF, with a 95% success rate in a single session [1, 19, 20].

Major Complications

Major complications include neurologic damage of contiguous structures and neurodystrophia (<0.5% of patients) [21].

Recurrence Rate

The recurrence rate for osteoid osteoma is approximately 5% (mainly in intra-articular lesions) [1].

RFA

Patient Selection and Technique

The indications for RFA in bone can be curative or palliative [22].

Major Indications

Curative indications include patients with primary (e.g., osteoid osteoma, osteoblastoma, chondroblastoma) or secondary bone tumours with contraindications to surgery or patient refusal of surgery (tumour size should be <5 cm in diameter). Palliative indications include painful bone metastases (pain not controlled by conventional therapies).

Ablation protocols vary according to lesion size and generator used. The best guidance modality is CAT with or without fluoroscopy [23]. Bone RFA requires regional block or general anaesthesia. At the end of ablation, 5 to 10 ml

ropivacaine 2 mg/ml is injected (strictly extravascularly) in contact with the periosteum to decrease postprocedural pain.

In osteoid osteoma, RFA technique is similar to that for laser ablation. The small size of the nidus does not require large ablation with expandable and multitined or perfusion and internally cooled electrodes. A 1-cm active-tip electrode is used to reach a temperature of 90°C within 6 to 10 minutes. If cortical bone must be passed through, a bone biopsy needle must be used, and the electrode is coaxially inserted into the lesion. Because the bone needle is not insulated, it should not come into contact with the active part of the RFA electrode. For tumours in close proximity to neurologic structures or other organs, thermal protection techniques are required [24]. RFA seems to be promising for the treatment of chondroblastomas, although few published data exist on outcomes exist [25–27].

For large volumes of ablation in tumours involving weight-bearing bones, additional consolidation with cementoplasty or surgery should be considered to prevent the risk of secondary fracture [28]. If cementoplasty is performed in the same session with RFA, the injection should be delayed until the tumour temperature has decreased to normal levels to avoid setting the cement too quickly.

The major advantage of RFA compared with chemical ablation (ethanol) is better delimitation of the ablation zone without risk of leak. The difficulty in treating bone tumours with RFA lies in thermal protection of vulnerable surrounding structures (particularly nerve roots) and penetration of the lesion, which in some cases requires coaxial introducer systems [2, 29].

Patient Outcome

For osteoid osteoma, the success rate of bone RFA is similar to that of laser therapy (>85% of patients), and the recurrence rate has been reported to be 5–10% [20, 30]. A 10–15% recurrence rate has been noted for typical osteoblastoma, and a 50% rate has been reported for aggressive osteoblastoma. An acute postablation inflammatory reaction can occur and should be prevented with administration of anti-inflammatory therapy. When using RFA for pain palliation in painful malignant bone tumours, significant (>50%) and rapid pain relief, including a subsequent substantial decrease in medication, has been reported in 70–95% of patients [4, 21, 31–34]. Because of advanced disease, recurrence of pain in relation to other metastases is common, but most patients remain pain free at the ablated area.

Complications

Complications include thermal damage of contiguous structures (insulation techniques and thermal monitoring are mandatory to decrease these risks) and infection.

Cryoablation

Patient Selection and Technique

The indications and procedural technique are similar to those for RFA. Cryoablation is efficient for both lytic and sclerotic tumours. Procedure planning remains the same. However, compared with RFA, ≤ 25 cryoprobes (17-gauge; best gap between probes is 2 cm) can be inserted simultaneously [35].

The procedure is performed with the patient under sedation or general anaesthesia using CAT or MRI guidance (MRI-compatible cryoprobes exist on the market). Percutaneous cryoablation appears to require less analgesia than RFA [36]. Cryoablation of skeletal metastases is a time-consuming procedure [37] because two 10-minute freeze cycles, separated by an 8-min passive thaw, should be performed per position. Monitoring of the ice ball (of predictable geometry based on the length and diameter of the expansion room at the tip of the probe) is necessary and can be achieved with CAT imaging (e.g., hypodense ice ball). For complete necrosis of the tumour, it is important to extend the margins of the ice ball by a minimum 5-mm distance beyond the tumour margins to ensure complete cell death.

When thermal insulation is required, fluid should be avoided. For the best tissue insulation and protection of surrounding organs, carbon dioxide (CO₂) should be used [24].

For large volumes of ablation in tumours involving weight-bearing bones, additional consolidation with cementoplasty or surgery should be considered to prevent the risk of secondary fracture [28, 38]. If cementoplasty is performed in the same session with cryoablation, injection should be delayed until the temperature of the tumour has increased to normal levels to avoid cement leakage.

Compared with other ablation techniques, cryotherapy offers the major advantages of direct visualization of the ice ball, simultaneous use of ≤ 25 probes, and less post-procedural pain [37, 39].

Effective cryoablation tissue injury depends on four criteria:

1. excellent monitoring of the process;
2. fast cooling to a lethal temperature;
3. slow thawing; and
4. repetition of the freeze–thaw cycle [35, 40–44].

Patient Outcome

Percutaneous cryoablation is a safe and effective method for palliation of pain caused by metastatic disease and has been reported to provide significant pain relief in $\leq 75\%$ of patients [37].

Complications

Complications include cryoshock phenomenon, a rare syndrome that happens after large ablations, especially in the liver, occurs in 1% to 3.8% of patients [45–47] as well as thermal damage to contiguous neural structures (temporary neuropraxia at -20°C and permanent neurologic damage at $\leq 40^{\circ}\text{C}$) [48].

Radiofrequency Ionization

Patient Selection and Technique

The best candidates for this technique are (1) patients having painful nonsurgical spinal tumours with intracanalicular extension and (2) patients with rupture of the posterior wall of the vertebral body and thus a high risk of cement leakage or shift of the tumour in the canal, both of which have neurologic consequences.

A bone trocar is used to allow tumoral access and coaxial insertion of the RFA 16-gauge bipolar electrode. A side-arm catheter is connected to the electrode to slowly inject saline solution for activation of the plasma field. The angled tip of the electrode allows the digging of several channels inside the tumour by rotation of the electrode on its axis. In spinal lesions after cavitation, acrylic cement is injected into the cavity for consolidation [21, 48].

Patient Outcome


Decreased neuralgia has been reported in $\leq 87\%$ of patients [48] (due to tumour decompression). Decreased cement leakage, although reported, requires further evaluation.

Complications

Complications include accidental damage to nerve roots and perforation of adjacent organs.

Thermal Insulation and Monitoring

Thermal insulation and monitoring is used to protect important anatomic structures from thermal injury during application of different thermal ablation techniques. Heating at $\geq 45^{\circ}\text{C}$ has been shown to be neurotoxic to the spinal cord and the peripheral nerves. Interposition of bone increases the insulation but depends on the thickness of bone lamella [2]. Precautionary measures should be taken when no intact insulating cortex remains between the tumour and the neurologic structures. The risk of thermal damage increases with the size of the active tip of the

Algorithm 1 Tumour-management therapeutic option algorithm

electrode [29, 49]. Thermal insulation techniques and thermal monitoring can be used alone or in combination, if necessary, to decrease unintended thermal injury to non-target structures.

Thermal Insulation

Thermal insulation consists of fluid, gas, or balloon interposition between adjacent nontarget structures and the ablation zone [1, 50–52]. Hydrodissection is one of the most commonly used techniques. With aerodissection, CO₂ is injected using a dedicated CO₂ injection syringe, which allows precise control of the volume of gas insufflated in the desired zone under positive pressure [24].

- *RFA* Saline is not a suitable solution (because of its high electrical conductivity), and dextrose in water 5% has been advocated.
- *Cryoablation* Carbon dioxide can be used to displace the nontarget structures away from the treatment zone, thus creating an excellent thermal insulator blanket. Fluid is not suitable because it can freeze in contact with the ice ball. In ultrasound-guided interventions, CO₂ dissection cannot be used (the ring-down artifact from gas deteriorates the image).

Thermal Monitoring

Thermal monitoring is achieved with thermocouples inserted coaxially through an 18-gauge spinal needle [53]. Special care is needed with RFA because the tip of metal-based thermocouple should never be placed too close to the RFA electrode (this could result in arcing and electrical conductivity).

Conclusion: Proposed Treatment Strategy

- Precise clinical evaluation of the patient is mandatory and should include origin and location of pain, previous treatment, patient-tolerated anaesthesia, and life expectancy.
 - Whole-body 3-D imaging should be used to precisely analyze the lesions and their relation to surrounding structures.
 - Multidisciplinary decision making is required to choose the most efficient and least disabling technique.
1. For a painful tumour involving a flat weight-bearing bone without invasion of surrounding tissues and associated with risk of compression fracture, single cementoplasty is the most appropriate technique.
 2. For a painful bone tumour with extension into surrounding soft tissues, thermal ablation is required to control the pain caused by soft-tissue invasion. However, if there is a risk of pathologic fracture, an additional consolidation technique is required (cementoplasty for flat bones or surgery for long bones).
 3. For nonsurgical spinal tumours extending toward the canal with rupture of the posterior wall and epidural extension, percutaneous tumour decompression using radiofrequency ionization is the best technique. After tumour decompression, the cavity created can be filled with cement (Algorithm 1).

References

1. Gangi A, Alizadeh H, Wong L, Buy X, Dietemann JL, Roy C (2007) Osteoid osteoma: percutaneous laser ablation and follow-up in 114 patients. *Radiology* 242:293–301

2. Dupuy DE, Hong R, Oliver B, Goldberg SN (2000) Radiofrequency ablation of spinal tumors: temperature distribution in the spinal canal. *AJR* 175:1263–1266
3. Goldberg SN (2001) Image-guided radiofrequency tumour ablation: challenges and opportunities. II. *J Vasc Interv Radiol* 12:1135–1148
4. Thanos L, Mylona S, Galani P et al (2008) Radiofrequency ablation of osseous metastases for the palliation of pain. *Skeletal Radiol* 37:189–194
5. Gangi A, Dietemann JL, Schultz A, Mortazavi R, Jeung MY, Roy C (1996) Interventional radiologic procedures with CT guidance in cancer pain management. *Radiographics* 16:1289–1386
6. Gangi A, Guth S, Imbert JP, Marin H, Dietemann JL (2003) Percutaneous vertebroplasty: indications, technique, and results. *Radiographics* 23:e10
7. Jang JS, Lee SH (2005) Efficacy of percutaneous vertebroplasty combined with radiotherapy in osteolytic metastatic spinal tumors. *J Neurosurg* 2:243–248
8. Weill A, Chiras J, Simon JM, Rose M, Sola-Martinez T, Enkaoua E (1996) Spinal metastases: indications for and results of percutaneous injections of acrylic surgical cement. *Radiology* 199:241–247
9. Gangi A, Kastler B, Klinkert A, Dietemann JL (1994) Injection of alcohol into bone metastases under CT guidance. *J Comput Assist Tomogr* 18:932–935
10. Baust J, Chang Z (1995) Underlying mechanisms of damage and new concepts in cryosurgical instrumentation. International Institute of Refrigeration, Paris
11. Mazur P (1997) The role of intracellular freezing in the death of cells cooled at supraoptimal rates. *Cryobiology* 14:251–272
12. Mathis JM WW (2003) Percutaneous vertebroplasty: technical considerations. *J Vasc Interv Radiol* 14:953–960
13. McGraw JK, Cardella J, Barr JD et al (2003) Society of Interventional Radiology quality improvement guidelines for percutaneous vertebroplasty. *J Vasc Interv Radiol* 14(Suppl):S311–S315
14. Yamada K, Matsumoto Y, Kita M, Yamamoto K, Kobayashi TTT (2004) Long-term pain-relief effects for patients undergoing percutaneous vertebroplasty for metastatic vertebral tumor. *J Anst* 18:292–295
15. Cotten A, Dewatre F, Cortet B et al (1996) Percutaneous vertebroplasty for osteolytic metastases and myeloma: effects of the percentage of lesion filling and the leakage of methyl methacrylate at clinical follow-up. *Radiology* 200:525–530
16. Shimony JS, Gilula LA, Zeller AJ, Brown DB (2004) Percutaneous vertebroplasty for malignant compression fractures with epidural involvement. *Radiology* 232:846–853
17. Buy X, Basile A, Bierry G, Cupelli J, Gangi A (2006) Saline-infused bipolar radiofrequency ablation of high-risk spinal and paraspinal neoplasms. *AJR Am J Roentgenol* 186(Suppl):S322–S326
18. Laredo JD, Hamze B (2004) Complications of percutaneous vertebroplasty and their prevention. *Skeletal Radiol* 33:493–505
19. Witt JD, Hall-Craggs MA, Ripley P, Cobb JP, Bown SG (2000) Interstitial laser photocoagulation for the treatment of osteoid osteoma. *J Bone Joint Surg Br* 82:1125–1128
20. Rosenthal DI, Hornicek FJ, Torriani M, Gebhardt MC, Mankin HJ (2003) Osteoid osteoma: percutaneous treatment with radiofrequency energy. *Radiology* 229:171–175
21. Gangi A, Guth S, Guermazi A (2009) Imaging in percutaneous musculoskeletal interventions. Springer, Berlin-Heidelberg, Germany
22. Gangi A, Basile A, Buy X, Alizadeh H, Sauer B, Bierry G (2005) Radiofrequency and laser ablation of spinal lesions. *Semin Ultrasound CT MR* 26:89–97
23. Gangi A, Kastler BA, Dietemann JL (1994) Percutaneous vertebroplasty guided by a combination of CT and fluoroscopy. *AJNR Am J Neuroradiol* 15:83–86
24. Buy X, Tok CH, Szwarc D, Bierry G, Gangi A (2009) Thermal protection during percutaneous thermal ablation procedures: interest of carbon dioxide dissection and temperature monitoring. *Cardiovasc Intervent Radiol* 32:529–534
25. Christie-Large M, Evans N, Davies AM, James S (2008) Radiofrequency ablation of chondroblastoma: procedure, technique and MR imaging follow up of four cases. *Skeletal Radiol* 37:1011–1017
26. Erickson JK, Rosenthal DI, Zaleske DJ, Gebhardt M, Cates J (2001) Primary treatment of chondroblastoma with percutaneous radiofrequency heat ablation: report of three cases. *Radiology* 221:463–468
27. Rybak DL, Rosenthal ID, James WC (2009) Chondroblastoma: radiofrequency ablation—alternative to surgical resection in selected cases. *Radiology* 251:599–604
28. Toyota N, Naito A, Kakizawa H et al (2005) Radiofrequency ablation therapy combined with cementoplasty for painful bone metastases: initial experience. *Cardiovasc Intervent Radiol* 28:578–583
29. Bitsch RG, Rupp R, Bernd L, Ludwig K (2006) Osteoid osteoma in an ex vivo animal model: temperature changes in surrounding soft tissue during CT-guided radiofrequency ablation. *Radiology* 238:107–112
30. Rosenthal DI, Springfield DS, Gebhardt MC, Rosenberg AE, Mankin HJ (1995) Osteoid osteoma: percutaneous radiofrequency ablation. *Radiology* 197:451–454
31. Callstrom MR, Charboneau JW, Goetz MP et al (2002) Painful metastases involving bone: feasibility of percutaneous CT- and US-guided radio-frequency ablation. *Radiology* 224:87–97
32. Kojima H, Tanigawa N, Kariya S, Komemushi A, Shomura YSS (2006) Clinical assessment of percutaneous radiofrequency ablation for painful metastatic bone tumors. *CVIR* 29:1022–1026
33. Callstrom MR, Charboneau JW (2007) Image-guided palliation of painful metastases using percutaneous ablation. *Tech Vasc Interv Radiol* 10:120–131
34. Goetz MP, Callstrom MR, Charboneau JW et al (2004) Percutaneous image-guided radiofrequency ablation of painful metastases involving bone: a multicenter study. *J Clin Oncol* 22:300–306
35. Theodorescu D (2004) Cancer cryotherapy: evolution and biology. *Rev Urol* 6(Suppl 4):S9–S19
36. Allaf ME, Varkarakis IM, Bhayani SB, Inagaki T, Kavoussi LR, Solomon SB (2005) Pain control requirements for percutaneous ablation of renal tumors: cryoablation versus radiofrequency ablation—Initial observations. *Radiology* 237:366–370
37. Callstrom MR, Atwell TD, Charboneau JW et al (2006) Painful metastases involving bone: percutaneous image-guided cryoablation—prospective trial interim analysis. *Radiology* 241:572–580
38. Bickels J, Kollender Y, Merimsky O, Isaakov J, Petyan-Brand R, Meller I (2004) Closed argon-based cryoablation of bone tumours. *J Bone Joint Surg Br* 86:714–718
39. Beland MD, Dupuy DE, Mayo-Smith WW (2005) Percutaneous cryoablation of symptomatic extraabdominal metastatic disease: preliminary results. *AJR Am J Roentgenol* 184:926–930
40. Baust JG, Gage AA (2004) Progress toward optimization of cryosurgery. *Tech Cancer Res Treat* 3:95–101
41. Clarke DM, Robilotto AT, Rhee E et al (2007) Cryoablation of renal cancer: variables involved in freezing-induced cell death. *Technol Cancer Res Treat* 6:69–79
42. Gage AA, Baust J (1998) Mechanisms of tissue injury in cryosurgery. *Cryobiology* 37:171–186

43. Hoffmann NE, Bischof JC (2002) The cryobiology of cryosurgical injury. *Urology* 60:40–49
44. Littrup PJ, Ahmed A, Aoun HD et al (2007) CT-guided percutaneous cryotherapy of renal masses. *J Vasc Interv Radiol* 18:383–392
45. Bageacu S, Kaczmarek D, Lacroix M, Dubois J, Forest J, Porcheron J (2007) Cryosurgery for resectable and unresectable hepatic metastases from colorectal cancer. *Eur J Surg Oncol* 33:590–596
46. Seifert JK, France MP, Zhao J et al (2002) Large volume hepatic freezing: association with significant release of the cytokines interleukin-6 and tumor necrosis factor α in a rat model. *World J Surg* 26:1333–1341
47. Seifert JK, Morris DL (1999) World survey on the complications of hepatic and prostate cryotherapy. *World J Surg* 23:109–113 (discussion 113–104)
48. Georgy BA, Wong W (2007) Plasma-mediated radiofrequency ablation assisted percutaneous cement injection for treating advanced malignant vertebral compression fractures. *AJNR Am J Neuroradiol* 28:700–705
49. Adachi A, Kaminou T, Ogawa T et al (2008) Heat distribution in the spinal canal during radiofrequency ablation for vertebral lesions: study in swine. *Radiology* 247:374–380
50. Diehn FE, Neeman Z, Hvizda JL, Wood BJ (2003) Remote thermometry to avoid complications in radiofrequency ablation. *J Vasc Interv Radiol* 14:1569–1576
51. Laeseke PF, Sampson LA, Brace CL, Winter TC III, Fine JP, Lee FT Jr (2006) Unintended thermal injuries from radiofrequency ablation: protection with 5% dextrose in water. *AJR Am J Roentgenol* 186(Suppl):S249–S254
52. Kam AW, Littrup PJ, Walther MM, Hvizda J, Wood BJ (2004) Thermal protection during percutaneous thermal ablation of renal cell carcinoma. *J Vasc Interv Radiol* 15:753–758
53. Nakatsuka A, Yamakado K, Takaki H et al (2009) Percutaneous radiofrequency ablation of painful spinal tumors adjacent to the spinal cord with real-time monitoring of spinal canal temperature: a prospective study. *CVIR* 32:70–75